


## STRUCTURAL CALCULATIONS

### Anchorage of Dental Equipment to Interior Wall: *CDX Backing over Steel Studs*

Prepared for: Henry Schein Dental



---

## TABLE OF CONTENTS

---

- CDX Attachment Analysis
- Lag Screw Analysis
- Backing Installation Detail
- Supporting Documentation

---

## SUMMARY OF RESULTS

---

- (2) 3/4" plies of CDX plywood shall be secured to existing cold-formed steel (CFS) stud framing using 2½" Simpson Strong-Tie #10 Self-Drilling Bugle-Head Screws (or approved equal) spaced at 12" maximum centers *per stud*. Specific Simpson Strong-Tie screw model number = F10T250BD.
- To adequately secure the KAVO Orthopantomograph OP 3D unit to wall/backing, use (2) 3/8" diameter x 2" long lag screws.
- To adequately secure the Planmeca Viso unit to wall/backing, use (4) 3/8" diameter x 2½" long lag screws.
- Refer to attached Backing Installation Detail for additional information.

# CDX ATTACHMENT ANALYSIS



CORRE, INC.  
Transportation & Municipal Engineering  
6510 Grand Teton Plaza, Ste 314  
Madison, WI 53719  
www.correinc.com

JOB HENRY SCHEN DENTAL

SHEET NO. CFS BACKING

OF

CALCULATED BY DETAILS ETP

DATE

1/4/22

CHECKED BY

DATE

SCALE NTS

### CALCULATE WITHDRAWAL CAPACITY OF CDX BACKING:

- RECOMMEND USING #10 SCREWS AT 12" CENTERS

→ PULL-OUT CAPACITY = 84 LBS/SCREW (SEE ATTACHED REFERENCE)


- CONSERVATIVELY ASSUME STUD SPACING = 24" [MAX SPACING PER 2018 INTERNATIONAL BUILDING CODE]

ASSUMES 8" WALL HT

$$\begin{aligned} \therefore \text{PULL-OUT CAPACITY} &= \left( \frac{9 \text{ SCREWS}}{\text{STUD}} \right) \left( \frac{48" \text{ BACKING WIDTH}}{24" \text{ STUD SPACING}} + 1 \right) \left( \frac{84 \text{ LBS}}{\text{SCREW}} \right) \\ &= 2,268 \text{ LBS} >> 1,400 \text{ LBS} \end{aligned}$$

OK ✓

### DETERMINE MINIMUM SCREW LENGTH:



### MIN SCREW LENGTH:

|                   |                         |
|-------------------|-------------------------|
| 1.500"            | CDX PLYES               |
| 0.120             | STUD THICKNESS          |
| 0.125             | SURFACE INCONSISTENCIES |
| 0.375             | MIN EXTENSION THRU STUD |
| $\Sigma = 2.120"$ |                         |

$\therefore$  SAY 2 1/2" LENGTH

Calculated withdrawal/pull-out capacity = 84 lbs/screw (Factor of Safety = 3.0)

REFERENCE

## Fastening Options

Conservatively assume lightest gauge steel studs

Recommend using #10 screws to secure (2) 3/4" CDX plies to the steel studs/framing.

Connections can be made using a variety of fastening options. It is critical to specify the proper fastener to ensure the proper performance of the connections in cold-formed steel construction. The most common and widely used connection methods are screw connections and weld connections. Each type of connection method has various advantages and disadvantages. Therefore, we provide data for both types so you can choose your preferred connection method.



**Self-Drilling Screws**—externally threaded fasteners with the ability to drill their own hole and form, or “tap,” their own internal threads without deforming their thread and without breaking during installation. These screws are high-strength, one-piece fasteners and are used if the connection of multiple thicknesses of 33mil steel or thicker. One of the more

common self-drilling screws is a #10-16 x 5/8 HWH SD, which indicates a #10 diameter shaft, 16 threads per inch, 5/8 length, hex washer head self-drilling screw.

**Fillet Welds**—used to make lap joints, corner joints and T-joint connections. As the illustration suggests, the fillet weld is roughly triangular in cross-section, although its shape is not always a right triangle or an isosceles triangle. Weld metal is deposited in a corner formed by the fit-up of the two members and penetrates and fuses with the base metal to form the joint.

**Flare Welds**—used to join rounded or curved pieces.

- A Flare Bevel groove weld is commonly used to join a rounded or curved piece to a flat piece.
- A Flare V groove weld is commonly used to join two rounded or curved parts.

## ALLOWABLE SCREW DESIGN VALUES (LBS/SCREW)

| Mils (Gauge) | Design thickness (in) | F <sub>y</sub> (ksi) | F <sub>u</sub> (ksi) | #8 Screw<br>0.164" Dia, 0.272" Head<br>Pss = 1000 lbs; Pts = 1545 lbs |           |          | #10 Screw<br>0.190" Dia, 0.340" Head<br>Pss = 1400 lbs; Pts = 1936 lbs |           |          | #12 Screw<br>0.216" Dia, 0.340" Head<br>Pss = 2000 lbs; Pts = 2778 lbs |           |          | 1/4" Screw<br>0.250" Dia, 0.409" Head<br>Pss = 2600 lbs; Pts = 4060 lbs |           |          |
|--------------|-----------------------|----------------------|----------------------|-----------------------------------------------------------------------|-----------|----------|------------------------------------------------------------------------|-----------|----------|------------------------------------------------------------------------|-----------|----------|-------------------------------------------------------------------------|-----------|----------|
|              |                       |                      |                      | Tension                                                               |           |          | Tension                                                                |           |          | Tension                                                                |           |          | Tension                                                                 |           |          |
|              |                       |                      |                      | Shear                                                                 | Pull-Over | Pull-Out | Shear                                                                  | Pull-Over | Pull-Out | Shear                                                                  | Pull-Over | Pull-Out | Shear                                                                   | Pull-Over | Pull-Out |
| 33 (20g)     | 0.0346                | 33                   | 45                   | 164                                                                   | 212       | 72       | 177                                                                    | 265       | 84       | 188                                                                    | 265       | 95       | 203                                                                     | 318       | 110      |
| 43 (18g)     | 0.0451                | 33                   | 45                   | 244                                                                   | 276       | 94       | 263                                                                    | 343       | 109      | 280                                                                    | 345       | 124      | 302                                                                     | 415       | 144      |
| 54 (16g)     | 0.0566                | 33                   | 45                   | 333                                                                   | 346       | 118      | 370                                                                    | 433       | 137      | 394                                                                    | 433       | 156      | 424                                                                     | 521       | 180      |
|              |                       | 50                   | 65                   | 333                                                                   | 500       | 171      | 467                                                                    | 625       | 198      | 569                                                                    | 625       | 225      | 613                                                                     | 752       | 261      |
| 68 (14g)     | 0.0713                | 33                   | 45                   | 333                                                                   | 436       | 149      | 467                                                                    | 545       | 173      | 557                                                                    | 545       | 196      | 600                                                                     | 656       | 227      |
|              |                       | 50                   | 65                   | 333                                                                   | 515       | 215      | 467                                                                    | 645       | 249      | 667                                                                    | 788       | 284      | 866                                                                     | 948       | 328      |
| 97 (12g)     | 0.1017                | 33                   | 45                   | 333                                                                   | 515       | 213      | 467                                                                    | 645       | 246      | 667                                                                    | 778       | 280      | 867                                                                     | 936       | 324      |
|              |                       | 50                   | 65                   | 333                                                                   | 515       | 307      | 467                                                                    | 645       | 356      | 667                                                                    | 926       | 405      | 867                                                                     | 1352      | 468      |
| 118 (10g)    | 0.1242                | 33                   | 45                   | 333                                                                   | 515       | 260      | 467                                                                    | 645       | 301      | 667                                                                    | 926       | 342      | 867                                                                     | 1143      | 396      |
|              |                       | 50                   | 65                   | 333                                                                   | 515       | 375      | 467                                                                    | 645       | 435      | 667                                                                    | 926       | 494      | 867                                                                     | 1353      | 572      |

### Screw Value Notes:

- 1 Capacities are based on Section J4 of the NASPEC AISI S100-16 Calculations.
- 2 When connecting materials of different steel thickness or tensile strengths, use the lowest values. Tabulated capacities assume two sheets of equal thickness are connected.
- 3 Calculated capacities are based on Allowable Strength Design (ASD) and include appropriate safety factors.
- 4 When multiple fasteners are used, screws are assumed to have a center-to-center spacing of at least 3 times the nominal diameter.
- 5 Screws are assumed to have a center-of-screw to edge-of-steel dimension of at least 1.5 times the nominal diameter of the screw.
- 6 Pullout capacity is based on the lesser of pull-out capacity in sheet closest to screw tip or tension strength of screw.
- 7 Pullover capacity is based on the lesser of pullover capacity for sheet closest to screw header, or tension strength of screw.
- 8 Listed capacities are for pure shear or tension loads only. For combined shear and pullover, see AISI Section J4.5.
- 9 Shear strength for #8, #10, #12, and 1/4" screws shall be greater than or equal to 1000 lbs, 1400 lbs, 2000 lbs, and 2600 lbs respectively.
- 10 Tension strength for #8, #10, #12, and 1/4" screws shall be greater than or equal to 1545 lbs, 1936 lbs, 2778 lbs, and 4060 lbs respectively.


Shear and Pull-Over loading conditions not applicable.

## ALLOWABLE WELD VALUES (LBS/IN)

| Mils (Gauge) | Design thickness (in) | F <sub>y</sub> (ksi) | F <sub>u</sub> (ksi) | Weld (1 inch) |              |
|--------------|-----------------------|----------------------|----------------------|---------------|--------------|
|              |                       |                      |                      | Fillet        | Flare Groove |
| 43 (18)      | 0.0451                | 33                   | 45                   | 619           | 544          |
|              |                       | 33                   | 45                   | 822           | 682          |
| 54 (16)      | 0.0566                | 50                   | 65                   | 1188          | 985          |
|              |                       | 33                   | 45                   | 1082          | 859          |
| 68 (14)      | 0.0713                | 50                   | 65                   | 1563          | 1241         |
|              |                       | 33                   | 45                   | 1269          | 1226         |
| 97 (12)      | 0.1017                | 50                   | 65                   | 1269          | 1402         |
|              |                       | 33                   | 45                   | 1550          | 1497         |
| 118 (10)     | 0.1242                | 50                   | 65                   | 1550          | 1712         |

### Weld Value Notes:

- 1 Capacities are based on Sections J2.5, and J2.6 of the NASPEC AISI S100-16 Calculations.
- 2 When connecting materials of different steel thickness or tensile strengths, use the lowest values. Tabulated capacities assume two sheets of equal thickness are connected.
- 3 Calculated capacities are based on Allowable Strength Design (ASD) and include appropriate safety factors.
- 4 Weld strengths are given in lbs/in and are based on E60XX electrodes.
- 5 For flare groove welds when t > 0.1", tw = (5/16)\*R when weld filled flush to surface.



Complies with AISI S100-16 • IBC 2018

# LAG SCREW ANALYSIS

Client: Henry Schein Dental

Project Description: CDX Backing Analysis (OP 3D unit)

Designed By: ETP

Date: 1/3/2022



# LAG SCREW WITHDRAWAL ANALYSIS


Analyzed in Accordance with the 2018 NDS Design Specifications (ASD & LRFD)

## LOADING INFORMATION:

|                                              |                       |          |                       |
|----------------------------------------------|-----------------------|----------|-----------------------|
| Withdrawal Load (ASD), P =                   | note: 1500N = 337 lbs | 340 lbs. | per Viso manual (p18) |
| Factored Withdrawal Load (LRFD), P =         | = 340 x 1.4           | 476 lbs. | Factored as noted     |
| Load Duration Factor, C <sub>D</sub> (ASD) = | Dead Load             | 0.90     | Table 2.3.2           |
| Time Effect Factor, λ (LRFD) =               | 1.4D                  | 0.60     | Appendix N; N.3.3     |

## LAG SCREW INFORMATION:

|                                               |                                |
|-----------------------------------------------|--------------------------------|
| Lag Screw Diameter, D =                       | 3/8"                           |
| Selected Length, L =                          | 2"                             |
| Root Diameter, D <sub>r</sub> =               | 0.265 in. Appendix A; Table L2 |
| Overall Thread Length, T =                    | 1.500 in. Appendix A; Table L2 |
| Unthreaded Length, S =                        | 0.500 in. Appendix A; Table L2 |
| Tapered Tip Length, E =                       | 0.219 in. Appendix A; Table L2 |
| Maximum Effective Thread Length = L - S - E = | 1.281 in.                      |
| Standard Washer Thickness =                   | 0.083 in.                      |
| No. of Lag Screws Provided =                  | 2                              |



## MEMBER GEOMETRY:

|                                                              |           |                                      |
|--------------------------------------------------------------|-----------|--------------------------------------|
| Side Member/Sheathing Thickness =                            | 1.000 in. | bracket + wall thickness             |
| Actual Thread Length Embedment Into Member, p <sub>t</sub> = | 0.698 in. |                                      |
| Specific Gravity of Member, G =                              | 0.50      | Average SG for plywood Table 12.3.3A |

## ANALYSIS:

Unadjusted Withdrawal Value, W = 1800 G<sup>3/2</sup> D<sup>3/4</sup> = 305 lbs/in. (12.2-1)

Adjusted Withdrawal Values:

For ASD: W' = W C<sub>D</sub> = 274 lbs/in.

For LRFD: W' = W 3.32 0.65 λ = 395 lbs/in.

No. Lag Screws Req'd:

ASD = P<sub>ASD</sub> / (p<sub>t</sub> x W'<sub>ASD</sub>) = 1.8 lags

LRFD = P<sub>LRFD</sub> / (p<sub>t</sub> x W'<sub>LRFD</sub>) = 1.7 lags

Check:

OK

OK

## RESULTS:

**Use a minimum of (2) 3/8" dia. x 2" lag screws to properly secure the OP 3D unit to the 2-ply 3/4" CDX plywood backing.**

Client: Henry Schein Dental

Project Description: CDX Backing Analysis (Viso unit)

Designed By: ETP

Date: 1/3/2022



# LAG SCREW WITHDRAWAL ANALYSIS


Analyzed in Accordance with the 2018 NDS Design Specifications (ASD & LRFD)

## LOADING INFORMATION:

|                                        |                      |            |                       |
|----------------------------------------|----------------------|------------|-----------------------|
| Withdrawal Load (ASD), $P =$           |                      | 1,400 lbs. | per Viso manual (p18) |
| Factored Withdrawal Load (LRFD), $P =$ | $= 1,400 \times 1.4$ | 1,960 lbs. | Factored as noted     |
| Load Duration Factor, $C_D$ (ASD) =    | Dead Load            | 0.90       | Table 2.3.2           |
| Time Effect Factor, $\lambda$ (LRFD) = | 1.4D                 | 0.60       | Appendix N; N.3.3     |

## LAG SCREW INFORMATION:

|                                                 |           |                               |
|-------------------------------------------------|-----------|-------------------------------|
| Lag Screw Diameter, $D =$                       | 3/8"      |                               |
| Selected Length, $L =$                          | 2-1/2"    |                               |
| Root Diameter, $D_r =$                          | 0.265 in. | Appendix A; Table L2          |
| Overall Thread Length, $T =$                    | 1.750 in. | Appendix A; Table L2          |
| Unthreaded Length, $S =$                        | 0.750 in. | Appendix A; Table L2          |
| Tapered Tip Length, $E =$                       | 0.219 in. | Appendix A; Table L2          |
| Maximum Effective Thread Length $= L - S - E =$ | 1.531 in. |                               |
| Standard Washer Thickness =                     | 0.083 in. |                               |
| No. of Lag Screws Provided =                    | 4         | <input type="text"/> override |



## MEMBER GEOMETRY:

|                                                     |           |                                         |
|-----------------------------------------------------|-----------|-----------------------------------------|
| Side Member/Sheathing Thickness =                   | 0.875 in. | bracket + wall thickness                |
| Actual Thread Length Embedment Into Member, $p_t =$ | 1.323 in. | CDX plies fully engaged                 |
| Specific Gravity of Member, $G =$                   | 0.50      | Average SG for plywood<br>Table 12.3.3A |

## ANALYSIS:

Unadjusted Withdrawal Value,  $W = 1800 G^{3/2} D^{3/4} =$  305 lbs/in. (12.2-1)

Adjusted Withdrawal Values:

For ASD:  $W' = W C_D =$  274 lbs/in.

For LRFD:  $W' = W 3.32 0.65 \lambda =$  395 lbs/in.

No. Lag Screws Req'd:

ASD  $= P_{ASD} / (p_t \times W'_{ASD}) =$  3.9 lags

LRFD  $= P_{LRFD} / (p_t \times W'_{LRFD}) =$  3.8 lags

Check:

OK

OK

## RESULTS:

**Use a minimum of (4) 3/8" dia. x 2-1/2" lag screws to properly secure the Viso unit to the 2-ply 3/4" CDX plywood backing.**

## BACKING INSTALLATION DETAIL



FOLLOW LOCAL CODES FOR ALL WORK.

DETAIL PROVIDED FOR BLOCKING & UTILITY REQUIREMENTS. REFER TO MANUFACTURER'S INSTALLATION GUIDE FOR ADDITIONAL INFORMATION.

MOUNTING MATERIALS SUPPLIED AND INSTALLED BY CONTRACTOR.

COLD-FORMED STEEL (CFS) WALL FRAMING. TYP.

CONFIRM HEIGHTS WILL NOT CONFLICT WITH X-RAY BRACKET HEIGHTS. SEE HENRY SCHEIN REP.



9C

PANORAMIC X-RAY

NOT TO SCALE

## SUPPORTING DOCUMENTATION

- Ensure the proper tightness and attachment of all of the screws and bolts when installing the device.
- All device covers must be properly installed before handing the device to the user.
- Ensure the electrical safety of the device by inspecting the proper grounding of user touchable metallic parts before handing the device to the user.

1500 newtons (N)  
= 337.2 lbs

## 2.5.2 Device installation requirements

### 2.5.2.1 Installation location requirements

**WARNING!** Ensure that each of the wall mount fixing screws and the wall can withstand pull-out force of at least 1500 N.


- The place where the device is to be installed and the position from where the user takes images must be correctly shielded from the radiation that is generated when the device is operated. Follow the local radiation and safety requirements.
- The device must be fixed to the wall and the floor.

**NOTICE!** If the device cannot be fixed to the floor, install the device to an exhibition stand. See chapter [Exhibition stand installation](#) on page 108 for instructions.

- The wall material should be suitable for fixing the device. If the wall is made of a weak material, you may have to use a reinforcing plate on the rear side of the wall to hold the fixing hardware.
- Make sure that the floor, where the device is to be installed, can support its weight. To avoid the device from tipping over, fix the device with floor bolts appropriate to the floor material. The floor bolts and the floor material must withstand pull-out force of at least 1500 N.
- The device must be installed on a hard floor surface. All soft and elastic material, such as a carpet, must be removed from under the column floor plate.
- Do not install the device in environments where corrosive or explosive vapours or flammable anaesthetics are present.
- Special steps regarding EMC need to be taken when installing the device. For more information, refer to the chapter [Electromagnetic Compatibility \(EMC\) tables](#) on page 149.
- The device is supplied with a 3 m (10 ft) long power cord. Ensure that the connected power cord and Ethernet cables are long enough, as they need to move along with the device's up/down movements.

**NOTICE!** It's recommended to route the power feed and Ethernet connection to behind the device to ensure unobstructed movement of the cables.

- It's recommended to use a max. 30 A circuit breaker with the device.
- Maximum allowed mains line impedance is 0.2 Ω.
- For permanent installation, a separate lockable mains switch (not supplied) is required to be installed to the mains feed.
- Recommended mains over-current releases:  
100-120 V: 16 A  
220-240 V: 10 A



1.b. Drill the holes to the wall.

Note the following requirements.

- If the wall is made of concrete or brick, use the M10x70 DIN 571 screws and the 14x70 expansion anchors to secure the wall bracket in position. Drill securing holes  $\varnothing 14$  mm (0.55 in.), 85 mm (3.3 in.) in depth, and insert the expansion anchors into the holes.
- If the wall is made of wood, use the M10x70 DIN 571 screws. Do not use expansion anchors with wooden wall. Drill securing holes  $\varnothing 7$  mm (0.3 in.), 70-75 mm (2.75-3 in.) in depth, for the mounting screws.

**NOTE**

Use four mounting screws (instead of two) if you use smaller anchors and/or screws than recommended.

**CAUTION**

The wall support installation must be able to resist a pull-out force of 6200 N (630 kg; 1400 lb).

1.c. Attach the wall support and tighten the screws firmly.